

Harmony™

PRENATAL TEST


An advanced blood test to assess the risk of fetal trisomies and evaluate the X and Y chromosomes

A simple, safe blood test

- * Highly accurate, individualized results for you, your practice and patients¹⁻⁶
- * Performed anytime after 10 weeks' gestation
- * Lowest cumulative false positive rate¹⁻⁶

Advanced Technology Behind the Harmony Test

The Advantages of Directed Analysis



cfDNA in blood

- Chr 21, 18, 13 cfDNA
- Other Chr cfDNA
- Unmapped cfDNA


- ▶ Cell-free DNA are short DNA fragments from chromosomes found in circulation
- ▶ In pregnancy, cell-free DNA from the fetus and mother are both present in maternal blood³
- ▶ The Harmony test uses efficient, directed analysis for accurate trisomy detection

Directed (Harmony Test)



Greater efficiency

MPSS (shotgun)



Random analysis of cfDNA


Informative Results

- ▶ Easy interpretation: Simple “High-Risk” or “Low-Risk” reporting for each trisomy
- ▶ Personalized risk results incorporate chromosome counts, fetal DNA fraction, gestational age, and maternal age
- ▶ 99.5% of risk score values are at either risk cutoff ($<1/10,000$ or $>99\%$) for autosomal trisomies
- ▶ X and Y Analysis with the Harmony Prenatal Test offers $>99\%$ accuracy for fetal sex⁷

Test Results

Fetal cfDNA Percentage: 10.5%


CHROMOSOME	RESULT	PROBABILITY	RECOMMENDATION
Trisomy 21 (T21)	HIGH RISK	Greater than 99/100 (99%)	Genetic counseling and additional testing
Trisomy 18 (T18)	Low Risk	Less than 1/10,000 (0.01%)	Review results with patient
Trisomy 13 (T13)	Low Risk	Less than 1/10,000 (0.01%)	Review results with patient
Y Analysis	Male Fetus	Greater than 99/100 (99%)	Review results with patient



Fetal fraction – key determinant for results

- ▶ A minimum amount of fetal cfDNA is necessary for reliable testing and quality results
- ▶ The Harmony test incorporates the measurement of fetal cfDNA into the analysis of **every sample**
- ▶ Increased maternal weight and early gestational age may contribute to the presence of low fetal cfDNA ($<4\%$)⁷

Effect of maternal weight as it correlates to $\geq 4\%$ fetal cfDNA (%)



Accurate:

Enhanced Performance with Individualized Results¹⁻⁷

Studied in over 6,000 patients, including >2,000 average-risk women¹⁻⁷

T21
T18
T13


Test Performance		False Positive Rate
T21	>99% (231 of 232)	<0.1%
T18	>98% (103 of 105)	<0.1%
T13	8 of 10	<0.1%

X and Y analysis is >99% accurate for fetal sex. It can also assess risk for sex chromosome conditions with test performance varying by the type of condition detected.⁷

- ▶ Only non-invasive prenatal test (NIPT) that has been exclusively evaluated in 1st trimester pregnant women
- ▶ Results in 99% of your patients with proper sample collection
- ▶ 95% of results reported within 9 days of accessioning⁷

Clinical Utility in a General Screening Population⁶

Nicolaides K.H., Syngelaki A., Ashoor G, et al., Noninvasive prenatal testing for fetal trisomies in a routinely screened first-trimester population. *Am J Obstet Gynecol* (2012); 207:374.e1-6.




- ▶ Both figures have the same number of patients
- ▶ 10 cases of trisomy 21 or trisomy 18
- ▶ 1,939 non-trisomy cases

Flexible for Multiple Patient Populations

- ▶ The Harmony Prenatal Test detects >99% of fetal trisomy 21 cases at a false positive rate of <0.1%
- ▶ Optional X and Y chromosome analysis available for fetal sex and X,Y sex chromosome analysis
- ▶ This test does not assess risk for mosaicism, partial trisomies or translocations
- ▶ The Harmony test is available for all singleton and twin pregnancies, including those conceived by IVF



Performance of Screening Tests for Trisomy 21^{4,7}



The Harmony Prenatal Test has been developed and is performed as a laboratory test service by Ariosa Diagnostics, a CLIA-certified clinical laboratory located in California, USA.

Ariosa™, Harmony™, and Harmony Prenatal Test™ are trademarks of Ariosa Diagnostics, Inc.
 ©2013 Ariosa Diagnostics, Inc. All rights reserved.
 Customer service: 1-855-9-ARIOSA (855-927-4672)

1. Sparks, A.B., Struble, C.A., Wang, E.T., Song, K., Oliphant, A., Non-invasive Prenatal Detection and Selective Analysis of Cell-free DNA Obtained from Maternal Blood: Evaluation for Trisomy 21 and Trisomy 18, *Am J Obstet Gynecol* (2012), doi: 10.1016/j.ajog.2012.01.030.
2. Ashoor, G., Syngelaki, A., Wagner, M., Birdir, C., Nicolaides, K.H., Chromosome-selective sequencing of maternal plasma cell-free DNA for first trimester detection of trisomy 21 and trisomy 18, *Am J Obstet Gynecol* (2012), doi: 10.1016/j.ajog.2012.01.029.
3. Sparks, A.B., Wang, E.T., Struble, C.A., Barrett, W., et al., Selective analysis of cell-free DNA in maternal blood for evaluation of fetal trisomy. *Prenat Diagn* (2012); 32(1):3-9. doi: 10.1002/pd.2922. Epub 2012 Jan 6.

4. Norton, M., Brar, H., Weiss, J., Karimi, A., et al., Non-Invasive Chromosomal Evaluation (NICE) Study: Results of a Multicenter, Prospective, Cohort Study for Detection of Fetal Trisomy 21 and Trisomy 18, *Am J Obstet Gynecol* (2012), doi:10.1016/j.ajog.2012.05.021.
5. Ashoor, G., Syngelaki, A., Nicolaides, K.H., et al., Trisomy 13 detection in the first trimester of pregnancy using a chromosome-selective cell-free DNA analysis method, *ULTRASOUND Obstet Gynecol* (2012), DOI: 10.1002/uog.12299.
6. Nicolaides K.H., Syngelaki A., Ashoor G, et al., Noninvasive prenatal testing for fetal trisomies in a routinely screened first-trimester population. *Am J Obstet Gynecol* (2012); 207:374.e1-6.
7. Internal data on file.