

Lung@MPD

| | | | |
|--------------|--------------------------------|---------------|---------------------------|
| <i>EGFR</i> | Hotspot mutation | <i>NTRK2</i> | Fusion |
| <i>ALK</i> | Fusion & Hotspot mutation | <i>NTRK3</i> | Fusion |
| <i>ROS1</i> | Fusion | <i>RET</i> | Fusion & Hotspot mutation |
| <i>KRAS</i> | Hotspot mutation | <i>FGFR1</i> | Fusion |
| <i>BRAF</i> | Fusion & Hotspot mutation | <i>FGFR2</i> | Fusion |
| <i>ERBB2</i> | Fusion, Hotspot mutation & CNV | <i>FGFR3</i> | Fusion |
| <i>MET</i> | Fusion, Hotspot mutation & CNV | <i>PIK3CA</i> | Hotspot mutation |
| <i>NTRK1</i> | Fusion | | |

Specimen:

Formalin-fixed paraffin embedded (FFPE) tissue

(Tumour content: >20% Tumour nuclei)

Turnaround time:

7-10 working days